



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG autoclave type VACUKLAV 31 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 2500 W,
Vacuklav 31 is an efficient steam autoclave offering 4 work programmes and 2 test programmes,
Product features:
Class B,
Automatic,
Pre-vacuum and fractional vacuum produced by vacuum pump,
Additional vacuum drying possible,
Preheat option helps maintain temperature between cycles,
Fitted with indicator lights to inform of unit errors,
2.5" LCD display,
Foiled control panel,
Keypad on the control panel in German,
Tray storage area on the autoclave lid,
Water drain connection on the front of the unit,
Sockets:
For Melaprint printer (printer not included),
For Meladem ion exchange system (system not included),
230 V socket,
Made of stainless materials,
Stable on legs, including front legs with height adjustment,
Has 4 work programmes and 2 test programmes:
Unpacked tools: pressure 2 bar - time 4 min - temperature 134 °C,
Packed tools: pressure 2 bar - time 5 min - temperature 134 °C,
Textiles: pressure 1 bar - time 17 min - temperature 121 °C,
Disinfection: pressure 0.3 bar - time 8 min - temperature 105 °C,
Bowie & Dick test: steam permeability test to check the autoclave's efficiency by the user,
Leakage test: vacuum test used to check the tightness of the autoclave,
Specification:
Sterilisation chamber capacity: 14,5 L,
Maximum load: 4 kg of instruments with trays or 1.5 kg of textiles,
Water consumption during sterilization: 300 - 400 mL, depending on the programme and loading level,
Included:
Screen trays with dimensions: 18 x 28 x 2 cm - 2 pieces,
Holder for trays,
Stand with holders for 3 trays,
Brand new condensate tank,
User manual in English (PDF),
Dimensions:
Autoclave: 63 x 45 x 45 cm,
Chamber: Ø 23 cm x 34 cm,
Weight: 42 kg,
It has a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN NEDERLANDS
МЫ ГОВОРИМ ПО РУССКИ
МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ
WE SHIP WORLDWIDE

