



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Kriopol R - Breeze II (30 l) Cryotherapy device (New)**

- New,
- Manufactured in Poland,
- Cryotherapy machine Kriopol R Bryza II is a simple and reliable device that allows effective and efficient cooling of the skin surface, provides full comfort and safety for the patient. Thanks to the manufacturer's many years of experience, we ensure the highest standard of the product,
- Advantages and benefits for users:
  - Pressure-free and safe to use,
  - Easy to use, convenient and reliable,
  - Selection of the tank depends on the needs of customers,
  - Unique technology allows the flexibility of the nitrogen vapor supply line to the patient's skin surface,
  - It is immediately ready for use - it reaches operating temperature (-100) to (- 160) degrees Celsius within 30 seconds of startup,
  - It does not cause skin frostbite, because the steam jet is uniform and does not contain liquid nitrogen droplets,
- Advantages and benefits for patients:
  - Reduction of pain,
  - Acceleration of metabolism,
  - Reducing the activity of the inflammatory process,
  - Reducing the use of anti-inflammatory and analgesic drugs,
  - Reduction of skeletal muscle tension,
  - Faster recovery and recuperation,
  - Improved mobility in the treated joints,
  - Acceleration of the healing process,
  - Reduction of swelling,
  - Reduces cellulite and firms the skin,
- Technical specifications:
  - Refrigerant: Liquid nitrogen (LN2),
  - Tank volume: 30 l,
  - Gas jet temperature: -160 °C (at the nozzle outlet),
  - Number of degrees of blowing intensity control: 4 + cryopuncture + 2 pulsation,
  - Liquid nitrogen consumption: 0.08 to 0.15 kg/min (continuous operation) depending on the set intensity of blowing,
  - Power supply: 230 V +/- 10 %,
  - Frequency: 50 Hz,
  - Protection class: I,
  - Fuse: 3.15 WTA,
  - Device dimensions (width / length / height): 470 / 540 / 1000 mm,
  - The device has an EC certificate in accordance with Directive 93/42/EEC Annex II and a declaration of conformity. It is notified to the register of medical devices,
- Valid technical passport issued,
- Warranty:
  - The device is covered by a 2-year manufacturer's warranty,
  - Dewar tanks - 12 months,
- Financing options (Poland only): Installments, Leasing, Loan

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

