



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEMMERT TV15U laboratory incubator (Reconditioned)

Used,
Technical condition: very good,
Visual condition: good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1500 W,
Memmert TV15U is a capacious and easy to use laboratory thermometer,

Equipped with:
Indicator lights for start-up and heater operation,
Knob for adjusting the temperature in the range: 50° - 220° C,
Built-in thermometer with a range of 0° - 250° C,
Two perforated trays with adjustable height for their position,

Dimensions:
Compartment: 25 x 47.5 x 26 cm,
Tray: 22.5 x 47 cm,
Overall: 43.5 x 59.5 x 56 cm,
Weight: 22 kg,

Comes with a current Technical Passport issued,

Warranty:
6 months for domestic market - Poland,
1 month for the international market,

Financing options (Poland only): Installments, Leasing, Loan,
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

