



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### Wooden chair for the waiting room (Reconditioned)

Used,  
Technical condition: good,  
Visual condition: good,  
Made in Poland,  
The chair is designed for waiting rooms in public utility rooms, e.g. in clinics, laboratories,  
It has visible signs of use,  
Easy to keep clean,  
Seats are set on a height of: 45 cm,  
Stable, on 4 legs with feet,  
Dimensions:  
Backrest: 43.5 x 39 cm,  
Seat: 44 x 37 cm,  
Overall: 50 x 49 x 80.5 cm,  
Weight: 4.5 kg,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for domestic market - Poland,  
1 month for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,  
If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

