



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Docking station (charger) AD ELEKTRONIK MBSS (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

German production,

Power supply: 230 V and 12 V,

Frequency: 50/60 Hz,

Power: 80 W,

AD Elektronik MBSS is a compact docking station for charging batteries (alternately),

It is designed to charge MEDTRONIC PHYSIO CONTROL batteries of the Fastpak, Fastpak2, Lifepak NiCd and Lifepak SLA types

Equipped with LEDs indicating charging, charging, conditioning and battery error,

Device panel - operation in three languages (German, English, French),

Included:

English manual (PDF),

Power cord,

Dimensions: 20 x 31 x 10 cm,

Weight: 3 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

