



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **PRAISTECH FIBI medical cabinet (grey) (New)**

New,  
Made in Poland,  
Fibi medicine cabinet has durable, metal construction,  
Equipped with double glazed doors with a lock,  
Four shelves provide ample storage space,  
Shelves can be adjusted in height thanks to perforated rails with holes every 1,5 cm,  
2 keys included,  
Dimensions:  
Shelves: 90 x 37.5 cm,  
General: 90 x 40 x 185 cm,  
Weight: 44.5 kg,  
Holds a current Technical Passport,  
Warranty: 24 months,  
Financing possibility (Poland only): Installments, Leasing, Loan

