



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Mobile trolley for MAQUET SERVO-U SERVO-n ventilator (Ex-Display)

New (demo),
Technical condition: perfect,
Visual condition: perfect,
Swedish production,
Powder coated aluminium alloy trolley for SERVO-U and SERVO-n series MAQUET ventilators,
Mobile on 4 double wheels 125mm GUY-RAYMOND with locking device,
Attachment base with locking device,
Equipped with transparent document compartment,
Dimensions: 64 x 55 x 55 cm,
Weight: 11.5 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan

