



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Microwave diathermy BOSCH RADARMED Type 12 S 250 (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 5.5 A,
Control panel in English,
Temporary work switch-off within the range: 0-30 minutes,
Infinitely adjustable DOSIS within: 0-250 (0-25),
F=2450 +/-50 MHz,
Applicator (electrode) with dimensions: 54 x 23 cm, mounted on a movable arm,
Mast height (together with electrode): 166 cm,
Mobile on four rubber wheels, two of them with brakes,
Overall dimensions: 66 x 80 x 148 cm,
Weight: 70,5 kg,
Cat 02,
Has a current Technical Passport issued,
Warranty:
6 months for the national market - Poland,
1 month for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan

