



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRAISTECH ECG trolley (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Made in Poland,

PRAISTECH trolley enables fast moving of the ECG apparatus, it is characterised by high durability and stability,

It has an upper table top and a lower basket which is removable and a comfortable handle for moving,

Height of the upper table: 85cm, lower basket: 50cm,

Mobile on four wheels with a diameter of 12cm, the front wheels swivel and the two rear wheels with a lock,

Dimensions:

General: 80 x 46 x 95cm,

Upper table top dimensions: 44 x 41 cm,

Lower basket measuring: 31 x 40 x 14 cm,

Weight: 25 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

