



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **HILL-ROM bedside cabinet with table (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
French production,  
Hill-Rom is a stylish bedside cabinet that can be used in medical centres, hospitals etc,  
Made of laminated chipboard,  
Ideal for rooms with limited space,  
Double sided opening drawers and doors allowing the cabinet to be placed on either side of the bed,  
Cabinet features:  
Easy to clean,  
Rounded edges on all sides of the cabinet, protect against damage,  
Features:  
Top with rounded edges measuring: 49 x 57 cm,  
Drawer equipped with pull-out separator, dimensions: 34.5 x 36.5 x 8.5 cm,  
Open shelf with dimensions: 42 x 42 x 17 cm,  
Lockable space divided into two levels - inside a sliding shelf, allowing the patient to adjust it to his/her needs,  
Foldable table with height adjustment from: 76 - 114 cm,  
A tilting top that allows the patient to position it at different angles for activities such as reading,  
The base of the cabinet is equipped with four double wheels,  
Dimensions:  
Table: 34 x 54 cm,  
Overall: 52 x 59/125 x 87 cm,  
Weight: 40 kg,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Installments, Leasing, Loan

