



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Endoscope trolley RICHARD WOLF 2191.331 (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Mobile endoscopic trolley Richard Wolf Model 2191.331 designed entirely for endoscopic equipment has two separate compartments and a lockable lower drawer,  
The upper top has two openings: for connecting cables ( $\varnothing$  6 cm) and for the monitor holder ( $\varnothing$  4 cm),  
Two Vertical ESG glass panels on the front are opened vertically with the possibility to hide them inside, they have magnetic closures,  
Inside the chambers, on the side walls, there are 6 IEC female power sockets (3 for each chamber),  
Possibility to lead the connection cables in the side wall column,  
Rear of trolley ventilated and removable,  
Two comfortable transport handles,  
Length of power cord: 470 cm,  
4 wheels with a diameter of 12 cm with a lock,  
Dimensions:  
General: 75 x 70 x 133 cm,  
Top: 62 x 57 cm,  
Shelves inside:  
Upper: 54 x 37 x 50 cm,  
Lower: 54 x 37 x 50 cm,  
Drawer: 51 x 36 x 12 cm,  
Weight: 125 kg,  
Comes with a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Installments, Leasing, Loan

