



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ADVAN AMM213TD medical monitor (Reconditioned)

Used,

Technical condition: very good,

Visual condition: very good,

Korean production,

Power supply: 230 V,

Frequency: 50/60 Hz,

The ADVAN AMM213TD is a professional medical monitor, which thanks to the variety of input and output ports can be used in many specialties in the field of medicine,

Specifications:

Screen: LCD,

Diagonal: 21.3 in,

Colour depth: 16.7 Million,

Aspect ratio: 4:3,

Response time:

Viewing angle: 85°,

Resolution: 1600 x 1200 @ 60Hz,

Input ports: DVI, D-SUB, Y-C, C-Video/SOG, SDI, R/Pr, G/Y, B/Pb, H/CS, VS,

Output sockets: Y-C, C-Video, SDI, R/Pr, G/Y, B/Pb, H/CS, VS

Menu language: English,

VESA mount with M4 x 15mm screws (stand and bracket not included),

Automatic detection of connected I/O devices,

Included:

English user manual (PDF),

24V power supply,

Power cable,

Dimensions: 51 x 39 x 9.5 cm,

Weight: 9 kg,

Comes with a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan

