



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**PHILIPS HD7 ultrasound scanner (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in U.S.A., (2012),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 780 VA,
- **Philips HD7** is an ultrasound system that combines advanced clinical features, excellent ergonomics and convenient digital data archiving tools,
- The device is equipped with:
  - HD7 Revision 3.1 system,
  - Monitor with easy to read 15" diagonal display,
  - English language keyboard on the control panel,
  - Menu selection in languages: EN, DE, FR, ES, IT, RU, PT, CN,
  - Integrated DVD station,
  - Linear head Model L12 - 03:
    - Frequency range: 3 - 12 MHz,
    - 15° trapezoidal imaging,
    - High resolution imaging for superficial applications including vascular, small area and musculoskeletal,
    - Supports 2D imaging, colour, Tissue Harmonic Imaging, PW Doppler and Color Power Angio,
- 6 head holders,
- Applications:
  - Abdominal,
  - Small organs and superficial,
  - Paediatric,
  - Musculoskeletal system,
  - Neonatal,
  - Urological,
  - Obstetrics,
  - Gynaecological and fertility,
  - Vascular,
  - Transcranial Doppler,
  - Cardiology,
- Imaging modes:
  - M - Mode,
  - CW (Continuous Wave) Doppler,
  - PW (Pulsed Wave) Doppler,
  - 2D,
  - Color Doppler,
  - CPA (Color Power Angio),
  - THI (Tissue Harmonic Imaging),
- Image archiving:
  - Ability to save images, image loops and test reports to the built-in hard drive and CD and DVD disks,
  - Windows compatible CD or DVD image storage formats,
  - Videoprinter (not included),
  - Possibility to directly connect a computer printer to the camera for printing examination reports,
- Equipped with the following inputs and outputs:
  - LAN,
  - S - Video,
  - S - VGA 9-pin,
  - D - Sub 15 - pin for footswitch connection,
  - 3 x USB 2.0 (including one at the control panel),
  - 2 composite video outputs,
- The device has a handy, convenient handles for movement,
- Mobile on 4 castors Ø 13 cm, 2 of them with lock and brake,
- Dimensions:
  - General: 98 x 55 x 136 cm,
  - Monitor: 41 x 33 x 7.5 cm,
- Weight: 77 kg,
- Cat 02,
- It has a valid Technical Passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

