



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS HD7 ultrasound scanner (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in U.S.A., (2012),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 780 VA,
- **Philips HD7** is an ultrasound system that combines advanced clinical features, excellent ergonomics and convenient digital data archiving tools,
- The device is equipped with:
 - HD7 Revision 3.1 system,
 - Monitor with easy to read 15" diagonal display,
 - English language keyboard on the control panel,
 - Menu selection in languages: EN, DE, FR, ES, IT, RU, PT, CN,
 - Integrated DVD station,
 - Linear head Model L12 - 03:
 - Frequency range: 3 - 12 MHz,
 - 15° trapezoidal imaging,
 - High resolution imaging for superficial applications including vascular, small area and musculoskeletal,
 - Supports 2D imaging, colour, Tissue Harmonic Imaging, PW Doppler and Color Power Angio,
- 6 head holders,
- Applications:
 - Abdominal,
 - Small organs and superficial,
 - Paediatric,
 - Musculoskeletal system,
 - Neonatal,
 - Urological,
 - Obstetrics,
 - Gynaecological and fertility,
 - Vascular,
 - Transcranial Doppler,
 - Cardiology,
- Imaging modes:
 - M - Mode,
 - CW (Continuous Wave) Doppler,
 - PW (Pulsed Wave) Doppler,
 - 2D,
 - Color Doppler,
 - CPA (Color Power Angio),
 - THI (Tissue Harmonic Imaging),
- Image archiving:
 - Ability to save images, image loops and test reports to the built-in hard drive and CD and DVD disks,
 - Windows compatible CD or DVD image storage formats,
 - Videoprinter (not included),
 - Possibility to directly connect a computer printer to the camera for printing examination reports,
- Equipped with the following inputs and outputs:
 - LAN,
 - S - Video,
 - S - VGA 9-pin,
 - D - Sub 15 - pin for footswitch connection,
 - 3 x USB 2.0 (including one at the control panel),
 - 2 composite video outputs,
- The device has a handy, convenient handles for movement,
- Mobile on 4 castors Ø 13 cm, 2 of them with lock and brake,
- Dimensions:
 - General: 98 x 55 x 136 cm,
 - Monitor: 41 x 33 x 7.5 cm,
- Weight: 77 kg,
- Cat 02,
- It has a valid Technical Passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

