



Water bath MEMMERT 200 VR (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 900 W,
- **Bath** made entirely of acid-resistant steel,
- Capacity: 4 L,
- Temperature control in the range: 20 - 110 °C,
- **The set includes:**
 - Drain hose,
- Dimensions:
 - External: 32 x 33 x 21 cm,
 - Internal: 20 x 20 x 11 cm,
- Weight: 6 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

