



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Laboratory incubator MEMMERT TV30U (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1500W,
- **The Memmert TV30U** is a capacious and easy-to-use laboratory incubator,
- Built-in thermometer with a range of 0 - 350° C,
- Smooth adjustment of temperature with a knob in the range: 30 - 300° C,
- Height-adjustable shelves,
- Fitted with two indicator lights showing when the unit is started and when the heater is operating,
- Adjustable fresh air inlet,
- **The set includes:** 2 shelves, dimensions 39.5 x 31 cm,
- Dimensions:
 - Overall: 51.5 x 50 x 68 cm,
 - Chamber: 37.5 x 37.5 x 33 cm,
- Weight: 25 kg,
- Cat 03,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

