



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**Anesthesia machine DATEX OHMEDA AESTIVA/5 (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in USA,
- Power supply: 230 V,
- Emergency power supply: via one new built-in battery: 12 V 2.3 Ah,
- Frequency: 50/60 Hz,
- AESTIVA 5 is an inhalation anesthesia machine used in operating, preoperative and recovery rooms,
- Has TWO vaporizer ports with Selectatec mounting system, no vaporizers included - can be retrofitted at additional cost,
- Has a ventilator including patient's circular cycle DATEX OHMEDA Model 7900 SMART VENT,
- Supply of gases: O2, N2O, Air,
- Selectable menu in language: PL, EN, ESP, FR, DE, IT, PT, NN, JA,
- Apparatus equipped with:
  - Compact patient breathing circuit, placed on a mobile arm, with oxygen sensor,
  - A suction machine connected to VAC,
  - Tabletop with dimensions: 31 x 47 cm,
  - Folding tabletop with dimensions: 32 x 27 cm,
  - Two drawers with dimensions: 39 x 27 x 14 cm,
  - Ergonomic handles for easy movement,
  - Space for two gas cylinders,
  - § Anesthetic gas scavenging system (AGS), hoses can be retrofitted at additional cost,
- The set includes:
  - GE Model D-FPD15-00 patient monitor/cardiac monitor with modular panel:
    - Used,
    - Technical condition: very good,
    - Visual condition: very good,
    - Taiwanese manufacture, (2010),
    - Power supply: 230 V,
    - Frequency: 50/60 Hz,
    - Power: 12 V,
    - Liquid crystal color LCD diagonal: 15",
    - Navigation knob makes it easy to enter information and navigate the screen,
    - Selectable menu in language: EN, DE,
    - Control panel in English,
    - Custom display options allow to view and analyze data in graphical and numerical formats, compile real-time measurements and trends, and freely configure each screen element - from curves to parameter names,
    - Audible and visual alarms,
    - Foiled control panel,
    - Equipped with the following modules:
      - GE Model E-CAi0V-00 anesthetic gas measurement with water trap:
      - GE Model E-PRESTN-00:
        - ECG - Brand new 3 - lead cable,
        - SPO2 - Brand new clip sensor,
        - NIBP - Brand new cuff with new air hose,
        - TEMP 1/2 - No sensor included (can be retrofitted at additional cost),
        - P1/P2 - No measurement line included (can be retrofitted at additional cost),
  - Gas supply hoses with customer's choice of tips (AGA or DIN),
  - Patient breathing circuit,
  - Sampling line for measuring gases,
  - Cardiac monitor user manual in English (PDF),
- Mobile, on four castors with central brake,
- Overall dimensions: 100 x 123 x 158 cm,
- Weight: 172.5 kg,
- In the domestic market (Poland), if the equipment is installed within 30 days of delivery, the installation and training are free of charge, after this deadline these services will be charged,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN DUTCH  
Мы говорим по русски  
Мы говорим по украински

