



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ERBE ERBOTOM T 130 Surgical diathermy (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 130 W,
- Infinitely adjustable power,
- Features the following functions:
  - Monopolar cutting - 130 W,
  - GRAD monopolar coagulation - 130 W,
  - Bipolar coagulation - 50 W,
  - Control panel in German,
- Included:
  - Neutral electrode cable,
  - Monopolar electrode holder with reusable cable for 4.0 mm diameter electrodes,
  - Bipolar cable, bipolar forceps can be retrofitted (at extra charge),
  - Brand new foot switch,
  - Power cord,
- Dimensions: 41 x 31 x 17 cm,
- Weight of set: 9.5 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

