



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ARJOHUNTLEIGH CARENDO Multifunctional hygiene chair (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photographs of the product,
- Made in Sweden,
- Power supply via 1 new built-in NiMH 24 V/2.5 Ah battery,
- Frequency: 50/60 Hz,
- Power: 150 VA,
- Carendo is a multifunctional hygiene chair that improves hygiene and showering routines, enabling the carer to work in a comfortable position,
- Designed so that a single carer can perform all hygiene tasks, including showering, toileting, washing hair, removing clothes, dressing, changing nappies and other tasks requiring discretion, with just a single patient transfer,
- It has an adjustable height and reclining backrest, so the carer and resident can set the best height and position for themselves during all hygiene procedures,
- The chair can be adjusted to three positions using the wired remote control:
  - Toilet seated (raising - high/low),
  - Semi-recumbent (shower),
  - Nursing (allowing access to the buttock area),
- Equipped with:
  - Neck cushion with adjustable height to suit patient height,
  - Full-length back rest, robust armrests and seat made of polypropylene,
  - Seat with toilet hole which allows overlap with toilet bowl,
  - Four swivel castors fitted with individual brakes,
  - Seat cushion (minor signs of use visible),
  - Emergency stop and lower seat,
  - Leg rest that automatically extends forward when the seat is adjusted to the semi-reclining position,
  - Audible low battery warning signal,
  - Comfortable handles to guide the chair,
  - Grid for storing e.g. toiletries at the back of the backrest,
- Included:
  - Batteries - 2 pcs,
  - External charger,
  - User manual in Polish (PDF),
- Dimensions:
  - Overall: 86.5 x 64 x 124 cm,
  - Seat height, lowest/highest position: 51 cm / 105 cm,
  - Lifting range (seat): 55 cm,
  - Min. length, lowest upright position: 87.5 cm,
  - Max. length, highest reclined position 151 cm,
  - Undercarriage: 86.5 x 64 cm,
  - Undercarriage, width between legs: 508 cm,
  - Minimum backrest height, lowest lying position 1020 mm
  - Turning diameter: Ø 96 cm,
- Maximum load capacity: 136 kg,
- Weight of the set: 64 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for the domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

