



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Hospital bed MMO Medical 4000 (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in France, (2020),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- MMO Model 4000 bed is designed for short-term care in medical institutions where medical care, supervision and medical control are provided when needed, and where electromedical devices used in medical procedures can be provided to maintain or improve the patient's condition. Such institutions include hospitals, nursing homes, rehabilitation centers, geriatric centers and other specialized care centers. Designed to allow comfortable sleeping and smooth entry and exit from the bed. This bed provides the best care for both patients and staff,
- Equipped with:
 - Two control panels for the patient with function locking and another for the caregiver with central locking,
 - Divided, folding side rails, easily lowered with a unique locking system,
 - Removable head and leg tops,
 - Two transport handles,
 - CPR handles on both sides,
 - 4 IV slots in each corner of the bed,
 - 4 corner bumpers,
 - Night light under the bed,
 - Integrated bed extension of 18 cm,
- Electric adjustment via control panel:
 - Backrest angle in the range: 0 - 70°,
 - Thigh support angle in the range: 0 - 30°,
 - Angle of calf support within the range: -20° to 20°,
 - Trendelenburg and anti - Trendelenburg position in the range: 14°,
 - Bed height in the range: 35 - 92 cm,
- CPR,
- Locking function,
- Mobile bed on 4 double wheels with a diameter of 15cm and a central brake,
- The set includes:
 - User manual in Polish (PDF),
- Dimensions: 214 x 100 cm,
- Maximum load capacity: 250 kg,
- Weight: 128 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

