



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Wooden waiting room bench with table (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Manufactured in Poland,
- The bench is designed for waiting rooms in public places, such as a clinic or laboratory,
- Equipped with:
  - Two seats with backs made of plywood, set at a height of: 46.5 cm,
  - Armrests at the ends of the bench,
  - Table with dimensions: 49 x 49 cm,
  - Stable base, on 4 chrome legs with feet,
- Dimensions:
  - Backrest: 43 x 38 cm.
  - Seat: 44 x 39 cm,
  - Overall: 58.5 x 174 x 83 cm,
- Weight: 24.5 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

