



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Laboratory Incubator ELKON (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Built-in digital controller BTC-9090:
 - Process value of display,
 - Setpoint display,
 - Alarm light,
- Maximum temperature: 95° C,
- Chamber equipped with:
 - 4 pcs. - Tray holders, with the possibility of hanging them at 14 height levels,
 - 2 pcs. - Screen trays with dimensions: 34,5 x 42 cm,
- Dimensions:
 - Chamber: 40 x 45 x 37 cm,
 - Overall: 63 x 72 x 55.5 cm,
- Weight: 65 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

