



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Strength Machine with weights stack TECHNOGYM KINESIS HIGH PULL (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Italy, (2014),
- TECHNOGYM Kinesis High Pull is designed to perform upper body pulling movements to strengthen muscles improving the efficiency of daily activities. It is designed with two adjustable roller pads that serve to provide the body with three-dimensional stability during pulling movements, while the footrests provide the user with additional stability while allowing the increased exercise variety,
- Specifications:
  - FullGravity technology is the only innovative technology that allows users to perform movements with resistance in all planes of motion,
  - Dual weight stack - each cable controls a single, independent level of resistance, allowing users to more effectively perform all movements, including alternating and reversing movements, based on the appropriate load for each side and training goal,
  - A patented cable system that rotates 360 degrees means users can perform all movements without the cables interfering with their forearms, unlike other traditional machines,
  - Kinesis allows to perform more than 400 exercises, free of charge and without restrictions, which can be performed completely alone or in the presence of a Personal Trainer,
  - The user can start exercising immediately, as the ergonomic handle is always ready to use and requires no prior settings or adjustments,
  - Innovative, attractive design,
  - Independently step adjustable load per side in 2.5 kg increments within the range: 2.5 - 35 kg,
  - Roller pads adjustable in 6 steps,
- It has normal signs of use including minor damage to the upholstery on the seat (the upholstery will be replaced with a new one, in the color chosen by the customer from our color chart),
- Overall dimensions: 128 x 165 x 187 cm,
- Maximum load: 180 kg,
- Weight: 320 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

