



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Biomak FE602 BIS Podiatry chair (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Poland,
- Power supply: 230 V,
- Power: 100W
- Podiatry chair BIOMAK FE 602 BIS provides comfort during the procedure not only for patients but also for those performing the procedure. Equipped with electric motors, it allows quick and convenient adjustment of the chair,
- Upholstery with medical certificate in salmon color (possibility to replace the upholstery with a new one from our color pattern for an additional fee),
- Specification:
 - 6 profiled cushions,
 - 5 actuators,
 - Controlled electronically by remote control,
 - Possibility of unfolding to recumbent position,
 - Infinite height adjustment in the range: 58 - 77 cm,
 - Tilt angle:
 - Backrest: 80°,
 - Seat: 20°,
 - Legs: 85°,
 - Ability to raise the back to an upright position,
 - Leg rest can be lowered,
 - Tilt of the entire seat,
 - Pilot holding handle,
 - Possibility of removing armrests,
 - Base on 4 legs,
- Includes:
 - Remote control,
 - Power cord,
- Dimensions:
 - Overall: 160 x 85 x 120 cm,
 - Length of seat unfolded to recumbent position: 185 - 200 cm,
- Weight: 55kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

