



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HILL-ROM AFFINITY TREE P3700 Delivery bed (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in USA,
- Power supply: 230 V,
- Emergency power supply: through two new built-in 12 V 3 Ah batteries,
- Frequency: 50/60 Hz,
- Hill-Rom Affinity Tree is a bed designed for use by women in labor,
- Electric bed adjustment via nurse and patient control panels:
 - Bed height with mattress in range: 56 - 97 cm,
 - Backrest adjustment from 0 - 63°,
 - Adjustment of the foot segment's height,
 - Inflating or deflating the lumbar mattress and seat for a certain firmness or softness,
- Manual adjustment:
 - Side-reversible, adjustable and removable footrests,
 - Removable foot section - removing transforms the bed into a gynecological chair,
 - Trendelenburg position: 8°,
- Equipped with:
 - Removable wooden headrest with built-in handles to help control and increase the mobility of the bed,
 - Folding headrest side rails with controls,
 - CPR lever - the quick-release segment lever is used especially in emergency situations (e.g. CPR), and allows quick leveling,
 - Folding delivery handles, covered with non-porous, fluid-resistant foam,
 - Air support system,
 - Adjustable knee pads,
 - Additional 12-volt AC accessory outlets (20 watts) are located on both sides of the bed,
 - Built-in night light - a photosensor allows automatic illumination, lighting up the center point of the floor and foot exit area,
 - Bowl for waste,
 - Bumpers,
- Mobile bed on casters with Ø 15 cm, with central brake,
- Includes:
 - Anti-decubitus mattress (has visible signs of use on the cover),
 - User manual in English (PDF),
 - Power cord,
- Dimensions:
 - Overall: 229 x 99 x 92.5/129.5 cm,
 - Delivery bed position (lying area): 198 cm long, 88.5 cm wide,
 - Gynecological chair position (lying area): 129 cm long, 87.5 cm wide,
- Maximum load capacity: 227 kg,
- Weight: 217 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

