



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HILL-ROM AVANTGUARD 1200 Hospital bed (ICU) (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in France,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 1.3 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 360 VA,
- AvantGuard 1200 is a multifunctional hospital bed designed to meet the needs of the entire medical team,
- Due to its features, it is designed for use not only in general care, but also in post-operative care and the ICU,
- It facilitates the use of vital function monitoring devices and patient transport,
- Has buttons for selective function locking located on both sides of the bed,
- Smooth electrical adjustment by means of a control panel on a flexible arm:
 - Bed height in the range: 55 - 95 cm,
 - Back segment tilt angle in the range: 0 - 65°,
 - Trendelenburg and anti - Trendelenburg position: + 16 / - 16°,
 - Thigh segment tilt angle: 36°,
 - Foot segment tilt angle: - 23°,
- Equipped with:
 - Safety folding handrails with holes that also serve as handles,
 - Back segment angle, Trendelenburg and anti-Trendelenburg positions indicator,
 - Central brake,
 - 4 drip spaces in each corner of the bed,
 - 4 corner bumpers,
 - Double-sided levers for CPR,
 - Side care unit with battery indicator (on both sides),
 - Accessory holder,
 - Removable head and leg tops,
 - 8 handles for attaching patient belts,
 - Integrated bed extension,
 - Brake lock alarm,
- X-RAY Transparency,
- Mobile, on 4 double wheels with Ø 15 cm, and an additional self-steering fifth wheel,
- The set includes:
 - User manual in English (PDF),
- Overall dimensions: 222 x 99 x 80/120 cm,
- Load capacity:
 - Safe: 185 kg,
 - Maximum: 250 kg,
- Weight: 175 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

