



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rotor for lower limb exercise (Reconditioned)

- Used,
- Technical condition: perfect,
- Visual condition: perfect,
- Rotor is designed for active and self-assisted exercises in the lower limbs, for people who suffer from muscle atrophy or contracture caused by long-term immobilization,
- It is also used in the rehabilitation of people who need to restore motor functions and re-learn the correct movement pattern lost through damage to the neurological system,
- The rotor has load adjustment by means of a knob, so it is possible to adjust it to the current functional status of limb and rehabilitation progress,
- The rotor can be used in hospitals, rehabilitation clinics and at home as a device for rehabilitation and individual recreation,
- Powder-coated steel frame,
- Equipped with:
 - Anti-slip caps,
 - Velcro fastening bands to secure the foot to the pedals,
- Dimensions: 48 x 34 x 40 cm,
- Weight: 3 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински



