



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BAUS AT Stretcher mount (New)

- New,
- Manufactured in Germany,
- Rail made of stainless steel,
- Equipped with:
 - 2 fixing points in the front part for attaching the front wheels of the stretcher,
 - 1 fixing point in the rear part with a lever for manual locking,
- Rail attached to the floor of the ambulance or to the base plate,
- Mount fits under the stretcher:
 - Stryker,
 - Ferno,
 - Stollenwerk,
 - Medirol,
- Overall dimensions:
 - Length: 180 cm,
 - Width: 26.5/38 cm,
 - Height 13.6 cm,
- Weight: approximately 11 kg,
- Valid technical passport issued,
- Warranty 24 months,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

