



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PENTAX ED-3480TK Video duodenoscope (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Japan,
- The ED-3480TK videoduodenoscope delivers - through a combination of Pentax DSP technology and high-resolution color CCD chips - endoscopic images of excellent quality in terms of resolution and color reproduction in a unique full-screen size,
- It is suitable for endoscopic examination of pediatric and adult patients,
- Specifications:
 - Working length of the viewing section: 1250 mm,
 - Minimum depth of field: 4 mm,
 - Outer diameter of the insertion tube: 11.6 mm,
 - Biopsy channel diameter: 4.2 mm,
 - Maximum tip bending angle:
 - Top: 120°,
 - Bottom: 90°,
 - Right: 110°,
 - Left: 90°,
 - Viewing angle: 80°,
- Compatible with:
 - Pentax EPK-1000,
 - Pentax EPK-700,
 - Pentax EPK- i,
- The picture shows minor damage to the endoscope, not interfering with its proper operation,
- Included:
 - Pentax case,
 - Valves - 2 pcs,
 - Lid,
- Dimensions of the case: 16 x 59 x 43 cm,
- Weight of the set: 4 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

