



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PENTAX EC-3890FK2 Video colonoscope (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Japan,
- Colonoscope EC-3890FK2 provide excellent high-resolution imaging in a full line of economical probes with ergonomic design and established quality,
- Includes high-resolution imaging sensors with bright illumination to ensure effective image clarity and clear details,
- Specifications:
  - Working length of the viewing section: 1500 mm,
  - Overall length: 1823 mm,
  - Depth of field: 3 - 100 mm,
  - Working part diameter: 13.2 mm,
  - Diameter of operating channel: 4.2 mm,
  - Maximum tip bending angle:
    - Top: 180°,
    - Bottom: 180°,
    - Right: 160°,
    - Left: 160°,
  - Viewing angle: 140°,
- Compatible with:
  - EPK i5000,
  - EPK-p,
- Includes:
  - Pentax case,
  - Valves - 2 pieces,
  - Lid,
- Dimensions of the case: 16 x 59 x 43 cm,
- Weight of the set: 4 kg,
- Cat 03,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

