



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PENTAX EC-3890Li Video colonoscope (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Japan,
- The EC-3890Li colonoscope is a reliable high-definition (HD) endoscope,
- It is equipped with high-definition (HD) CCD image sensors. The HD image quality ensures the clarity and sharpness of the endoscopic image, which can help the operator to visualize the gastrointestinal mucosa in detail,
- Technical specifications:
 - Working length of the viewing section: 1700 mm,
 - Overall length: 2023 mm,
 - Minimum depth of focus: 4 - 100 mm,
 - Insertion tube diameter: 13.2 mm,
 - Biopsy channel diameter: 3.8 mm,
 - Maximum tip bending angle:
 - Top: 180°,
 - Bottom: 180°,
 - Right: 160°,
 - Left: 160°,
 - Viewing angle: 140°,
- Includes:
 - Pentax case,
 - Valves - 2 pieces,
 - Lid,
- Dimensions of the case: 16 x 59 x 43 cm,
- Weight of the set: 4 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

