



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ELBANTON ET 50 Hypothermie Laboratory incubator (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Manufactured in Holland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- ELBANTON Model ET 50 is a large, capacious and easy to use laboratory thermometer,
- Temperature range: from ambient to + 50°C,
- Infinitely adjustable temperature with a knob,
- Built-in digital thermometer,
- Housing and chamber made of stainless steel,
- Possibility of manual adjustment of the drawers position height and the shelf on 13 levels (any arrangement),
- Adjustable drawer dividers (allows the operator to divide the drawer into 3 equal compartments),
- The set includes:
 - Mesh shelf,
 - Two mesh drawers,
 - Power cord,
- Dimensions:
 - General: 175 x 74 x 75 cm,
 - Compartments: 52 x 75 x 49 cm
 - Mesh shelf: 48 x 50 cm,
 - 2 x mesh drawer: 49 x 45 cm,
- Weight: 150 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

