



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DAUPHIN EDDY 04110.087 Cantilever chair (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Germany,
- Dauphin chair is designed for offices, cabinets, waiting rooms,
- It has visible minor traces of use,
- Plastic perforated backrest in gray color,
- Upholstered seat, set at a height of: 47 cm,
- Plastic armrests in black,
- Chrome cantilevered tubular frame,
- Dimensions:
 - Backrest: 46.5 x 33 cm,
 - Seat: 44 x 45.5 cm,
 - Overall: 60 x 58.5 x 80 cm,
- Weight: 9.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

