



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STIEGELMEYER 5085-04-85-00-0 Rehabilitation bed (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Stieglmeyer is a comfortable rehabilitation bed with electric height adjustment,
- Two-section,
- The bed can be retrofitted with a gofer mattress at an additional cost,
- Sturdy and stable construction ensures safety and comfort,
- Equipped with:
 - Wired remote control,
 - Wooden panels on head and leg sides with cutouts for convenient bed movement,
 - Folding handrails on both sides of the bed,
 - Bumpers in each corner,
- Electric, infinite adjustment by wired remote control:
 - Bed height in the range: 40 - 80 cm,
 - Angle of the back segment in the range: 68°,
- Bed and remote control have visible signs of use,
- Mobile, on 4 wheels with two central brakes,
- Dimensions:
 - Bed: 193 x 90 cm,
 - Overall: 209 x 102 x 90/130 cm,
- Weight: 126 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

