



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DORR Heating cabinet (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Holland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1300 W,
- DORR blanket heating cabinet is suitable for heating blankets and molton. Heating cabinets are used in areas where the patient can benefit from a little extra heat, such as wards for waking and nursing, but also dialysis rooms. Thanks to internal circulation and a perforated bottom, optimal convection takes place, so blankets and moltons are heated quickly and evenly. The chamber is equipped with a ventilation space covered by a stainless steel perforated sheet. As a result, proper convection also occurs when the heating chamber is fully loaded,
- Temperature control in the range: 34 - 40°C,
- Equipped with 5 shelves (height adjustable),
- Cabinet made of stainless steel,
- Base on 4 legs,
- The set includes:
  - Power cord,
  - 5 shelves,
- Dimensions:
  - Overall: 80 x 40 x 140 cm,
  - Chamber: 115 x 72 x 28 cm,
  - Shelf: 27 x 72 cm,
- Weight: 90 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

