



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BTL 1500 Gynecological chair (Ex-Display)

- New, post-lease,
- Technical condition: perfect,
- Visual condition: very good,
- Manufactured in Czech Republic,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 300 W,
- 3 - motor,
- BTL 1500 is used for gynecological and urological examinations and procedures,
- Upholstery with medical approval, orange color, has minor rubbing in the back segment (upholstery can be replaced with a new one, from our color pattern for an additional fee),
- Electric adjustment by wired remote control:
 - Angle of the back segment in the range: 10° - 55°,
 - Seat height in the range: 80 - 100 cm,
 - Seat tilt angle in the range: 0° / +45°,
- The chair has:
 - Headrest,
 - Patient handles,
 - Rails for attaching additional accessories,
 - Stirrups,
 - Retractable stainless steel metal bowl for waste, dimensions: 13.5 x 50.5 x 15 cm,
 - Hygienic paper holder for the seat,
 - Folding footrest with dimensions: 29.5 x 49.5 x 19 cm,
- The set includes:
 - Wired remote control,
- Dimensions:
 - Back segment: 90 x 60 x 6 cm,
 - Thigh segment: 42.5 x 61 x 6 cm,
 - Overall: 140 x 80 x 152 cm,
- Maximum load capacity: 180 kg,
- Weight: 80 kg,
- Cat 03,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

