



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KODAK MEDICAL X-RAY PROCESSOR 102 Processor (Reconditioned)

- Used, very good technical condition,
- Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50-60 Hz,
- Power: 2.2 kWh,
- Film transport: Continuous transport by rollers,
- Formats of processed films:
 - Minimum sheets 10 x 10cm,
 - Maximum sheet width 35.8cm,
- Processing capacity:
 - 24cm x 30cm: 89 films/hour,
 - 35cm x 43cm: 55 films/hour,
- Standard processing cycle is: 105 seconds,
- Adjustable developer temperature: 28 - 37°C,
- Water consumption: 1.9l/minute during the process,
- Anti-crystallization and anti-oxidation process,
- Noise level: 58dB(A),
- Weight: Empty 35kg, with solutions 50kg,
- Size: 77cm x 59cm x 42cm (112cm on the work table),
- User manual in Polish (PDF format),
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

