



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MIDMARK PROMOTAL ALIZEO 21765-10 Diagnostic and examination table (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in France,
- Power supply: through two new built-in batteries: 12 V 4.5,
- Frequency: 50/60 Hz,
- Power: 30 Watt,
- 2 - segment,
- Waterproof mattress with a thickness of 8 cm, medically approved, blue in color, has slight rubbing on the cover corners,
- Reinforced stretcher construction,
- Electric adjustment:
 - Bed height in the range: 64 - 98 cm,
- Manual adjustment:
 - Backrest angle by gas spring in the range: from 0 to 80°,
 - Trendelenburg 8° and ant-Trendelenburg 14° position adjustment,
 - Folding side rails,
- Equipped with:
 - Rotating corner bumpers,
 - Folding telescopic infusion rod (2 hooks) with height adjustment in the range: 62 - 129 cm,
 - Folding side rails,
 - Folding footrest handle, and handles on the end of the head,
 - Rails for attaching additional accessories,
 - Side basket with dimensions: 90 x 16.5 x 10 cm,
 - 5-The directional wheel,
 - Large wheels with a diameter of 20 cm, with a central brake,
- Translucent to X-Ray,
- Dimensions:
 - Bed: 205 x 64 cm,
 - Overall: 228 x 83.7 x 95/130 cm,
- Maximum load capacity: 200 kg,
- Weight: 129 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

