



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **POLYVAC Sterilization tray (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in USA,
- PolyVac sterilization tray is used to organize and store various accessories when not in use and to hold accessories during sterilization,
- Two clips protect the container from accidental opening,
- Dimensions: 17.5 x 29 x 5.5 cm,
- Weight: 0.6 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN DUTCH  
Мы говорим по русски  
Мы говорим по украински

