



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Contour 3D MEDTRONIC T7690 Accessory tray (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in USA,
- The Medtronic accessory tray is used to organize and store various accessories when not in use and to hold accessories during sterilization. The tray itself is not a sterile barrier. The accessories and tray must be placed in a sterilization pouch or autoclave sterilization cloth before sterilization,
- Two clips secure the container from accidental opening,
- Dimensions: 16 x 30 x 7.7 cm,
- Weight: 0.4 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

