



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **HETTICH UNIVERSAL 32 Laboratory centrifuge (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 400 VA,
- HETTICH UNIVERSAL 32 perfectly smooth operation at both high and low speeds, versatile and easily replaceable accessories, and compact size,
- In the low speed range (max. RCF 4193), it accommodates standard test tubes up to 100 ml, microtiter plates, Falcon test tubes, blood collection tubes and cyto accessories. Microlyte test tubes can be accelerated to max. RCF 23907. Centrifuge samplers with larger volumes (up to 85 ml) also benefit from high motor power. With the RCF 9418, bacteria and yeast sediment quickly and easily,
- Extremely short start-up and total stopping times, even with rotors that reach higher speeds,
- Wide range of accessories including pivoting rotors, angled rotors, cups and adapters for a wide variety of centrifuge tubes,
- Centrifugation chamber made of stainless steel,
- Automatic rotor recognition,
- Easily replaceable rotors,
- Emergency lid lock release,
- Motor overheating protection,
- Specifications:
  - RPM/RCF speed: 18000 / 23 907
  - Max. capacity: 4 x 100 ml,
  - Max. density: 1.2 kg/dm<sup>3</sup>,
  - Kinetic energy: 6500 Nm,
  - Ambient temperature: 2°C to 40°C
  - Relative humidity: max. 80% to 31°C, falling linearly to 50% at 40°C,
  - Noise level: 66 dB(A) depending on rotor,
  - Can be retrofitted with selected accessories (cups / adapters / samplers) at an additional cost (consider rotational speeds when selecting accessories),
- Included:
  - Hettich Rotor 1624,
  - User manual in English (PDF),
  - Power cord,
- Dimensions: 30 x 42 x 49 cm,
- Weight: 24 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

[illegible]