



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MASIMO LNCS NEO 2329 SPO2 sensor (New)

- New,
- Manufactured in USA,
- The LNCS sensor is designed for use with apparatuses:
 - Equipped with Masimo SET® oximetry functions,
 - Licensed for use with LNCS sensors,
 - Nellcor pulse oximeters or Nellcor-compatible pulse oximeters, except Nellcor OxiMax®,
- Specifications:
 - Neonatal, pediatric,
 - Sensor attachment method: adhesive,
 - Dedicated for single patient use,
 - Sterile,
 - Designed for:
 - up to 3 kg (adhesive attachment to the leg) from 40 kg (adhesive attachment to the finger),
 - Cable length 45 cm,
 - Latex-free,
 - Replaceable adhesive (patch),
 - Expiration date: 01-11-2024,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

