



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS GIF XQ-10 Gastrofibroscope (Reconditioned)

- Used,
- Technical condition: good - the gastrofibroscope is leaking,
- Visual condition: good,
- Manufactured in Japan,
- Viewing direction: forward view,
- The camera features no blind spots,
- Depth of field: 3-100 mm,
- Max. bend angle of the endoscope tip up and down: 190°/90°,
- Max. tip bending angle to the right and left: 100°/100°,
- Maximum field of view angle: 120°,
- Included:
 - Original case with dimensions: 59 x 43 x 16 cm,
 - Valves,
- Dimensions:
 - Outer diameter: 10.5 mm,
 - Outer diameter of the tip: 9.8 mm,
 - Working channel diameter: 2.8 mm,
 - Working length: 103 cm,
 - Total length: 133 cm,
- Weight of the set: 4.1 kg,
- Valid technical passport issued,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

