



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEMMERT U30 Laboratory incubator (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1400 W,
- MEMMERT U30 is a simple, compact and easy to use laboratory thermometer,
- Specifications:
 - Temperature control in the range: 30°C to 220°C,
 - 6-stage ventilation control,
 - Built-in thermometer,
 - Built-in timer,
 - Made of acid-resistant steel,
 - Chamber equipped with a mesh shelf,
- Dimensions:
 - Overall: 51 x 48 x 68 cm,
 - Chamber: 37 x 38 x 33 cm,
 - Shelf: 31 x 39 cm,
- Weight: 33 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

