



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**MEDIMA P300 Infusion (volumetric) pump (Reconditioned)**

- Used,
- Technical condition: perfect,
- Visual condition: perfect,
- Manufactured in Poland, (2020),
- Power supply: 230 V,
- Emergency power supply: through one built-in Ni - MH 9.6V / 2.0 Ah battery, (can be replaced with a new one at additional cost),
- Frequency: 50/60 Hz,
- Power: 20 VA,
- Medima P300 - the pump uses advanced solutions that results in safety. It is not only attractively priced, but also offers the best solutions and build quality,
- Product features:
  - Large 3.2" color graphic display with high resolution and contrast,
  - Intuitive user interface,
  - Touch screen, symbolic and alphanumeric keyboard,
  - Device menu in Polish,
  - Enhanced alarm system,
  - Automatic infusion set installation system,
  - Very short time to start infusion,
  - Exceptional infusion uniformity,
  - Multiple infusion modes, including continuous, intermittent, profile (up to 24 cycles),
  - Advanced and intuitive to use drug library or label library,
  - Wide selection of infusion types and units (ng, µg, mg, µU, mU, U, kU, µmol, mmol, mol, / kg / min, hr, 24 hr),
  - Advanced solutions resulting in infusion safety,
  - Automatic and manual bolus (ability to program dose, time or rate of delivery),
  - Automatic elimination of occlusion bolus,
  - 12 levels of occlusion pressure,
  - Drug library (120 drug dosing procedures with all necessary infusion parameters can be stored in the pump's memory, additionally, limitations can be added to each parameter, which warn about exceeding permissible values during infusion programming, with soft limits warning about exceeding recommended values, and hard limits prevent the entry of values outside their range, once the drug is selected from the list, the pump is ready for infusion, the drug library can be entered directly from the pump keyboard, transferred from the pump to the pump via RS232 interface, created in Medima Drug Editor software and transferred to the pump directly from the PC or via MedimaNet),
    - Integration with hospital databases (HIS/PDMS),
- Built-in outlets:
  - 12 V,
  - RS - 232,
  - To call personnel,
- Specifications:
  - Flow rate:
    - 0.1 - 99.99 ml/h with increments of 0.01 ml/h,
    - 100 - 999.9 ml/h with increments of 0.1 ml/h,
    - 1000 - 1200 ml/h with increments of 1 ml/h,
  - Volume and infusion time:
    - 0.1 - 20000 ml; 1 min - 200 hours,
  - Bolus, initial bolus:
    - Rate: 20 - 1200 ml/h,
    - Volume: 1 - 500 ml,
    - Automatic and manual bolus,
  - Occlusion:
    - 12 levels (75 ÷ 900 mmHg),
- Included:
  - Power cord,
- Overall dimensions: 20 x 28, 7 x 11.5 cm,
- Weight: 2.5 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.  
WE SPEAK ENGLISH  
HABLAMOS ESPANOL  
WIR SPRECHEN DEUTSCH  
ON PARLE FRANCAIS  
WIJ SPREKEN DUTCH  
Мы говорим по русски  
Мы говорим по украински

