



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PENTAX Endoscope cart (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in England,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Equipped with:
 - Top,
 - 4 shelves,
 - Power strip with 5 power outlets,
 - Side panels for power cords,
 - Rubber bumpers,
 - Convenient handles for moving,
 - 1000 VA isolation transformer,
 - Made of powder-coated steel,
 - Openable rear enclosure with key lock,
 - Mobile on 4 anti-electrostatic casters with a diameter of 15 cm including 2 with a brake,
- Includes:
 - Power cords (5 pieces),
- Dimensions:
 - Overall: 144 x 75 x 62 cm,
 - Top table: 49 x 50 cm,
 - Top shelf: 35 x 48 cm,
 - Middle shelves: 48 x 50 cm,
 - Bottom shelf: 41 x 48 cm,
- Weight: 94 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

