



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **SCHMITZ Drip stand (Reconditioned)**

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in Germany,
- Designed to hang containers of intravenous fluids for patient administration,
- Stand made of chrome-plated steel tubes,
- Telescopic design for easy manual height adjustment in the range: 120 - 213 cm,
- Equipped with a hanger featuring 4 hooks and a handle with a condensation container,
- Mobile on 5 - pointed mobile base with double casters including two with brake,
- Dimensions:
  - Maximum height: 215 cm,
  - Base diameter: 61 cm,
- Weight: 4 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

