



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NITROCARE NTCR SD07 Dialysis chair (Ex-Display)

- New,
- Technical condition: perfect,
- Visual condition: perfect,
- Manufactured in Turkey, (2021),
- Power supply: 230 V,
- Emergency power supply: through two new built-in batteries: 12 V 1.3 Ah,
- NITROCARE chair is designed for procedures such as dialysis and blood collection,
- 4 - segmented,
- Height adjustment by wired remote control:
 - Chair height by hydraulic actuator in the range: 65 - 105 cm,
 - Backrest angle: 72°,
 - Angle of leg segment: 36°,
- Manual tilt capability:
 - Trendelenburg position within 13°,
 - Anti-Trendelenburg position within 12°,
- Height and angle of armrests,
- Equipped with:
 - Manually height-adjustable serum hanger with single hook,
 - Head pillow,
 - Foot support to prevent patient moving,
 - Folding lockable 2 side armrests,
 - 4 bumper discs,
 - Drive system with four wheels with a diameter of 14.5 cm, and double-sided central locking mechanism,
- Dimensions:
 - Bed:
 - Head cushion: 23.5 x 5.5 x 27 cm,
 - Back segment: 60 x 7 x 67 cm,
 - Thigh segment: 66 x 7 x 46.5 cm,
 - Leg segment: 60 x 7 x 72 cm,
 - Foot support: 43.5 x 4.5 x 23 cm,
 - Overall: 218 x 65/105 x 90 cm,
- Maximum load capacity: 175 kg,
- Weight: 106.5 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

