



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LOVIBOND TC 135 S Laboratory incubator (Ex-Display)**

- New, ex-display,
- Technical condition: perfect,
- Visual condition: very good,
- Manufactured in Germany, (2021),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 270 VA,
- Lovibond TC 135 S is a thermostatic cabinet for temperature-controlled sample storage and BOD determination in wastewater analysis,
- Highly efficient components and reinforced insulation ensure low power consumption values,
- Temperature control between 2° C and 40° C, in steps of 0.1° C,
- Equipped with:
  - Illuminated LED display with set and current temperatures,
  - Convenient internal 230 V power sockets allow the use of mixing units or other electrical appliances without opening the door and lowering the temperature,
  - Built-in axial fans: 320 m3 /h capacity,
  - Door with built-in key lock,
  - 3 shelves with adjustable placement height,
  - Capacity: 135 L,
- Includes:
  - Key,
  - User manual in English (PDF),
  - Abbreviated user manual in Polish,
- Dimensions:
  - Chamber: 44.1 x 51.3 x 70.2 cm,
  - Overall: 60 x 60 x 85 cm,
- Weight: 39 kg,
- Valid technical passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

