



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HILL-ROM 900 Hospital bed (ICU) (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: very good,
- Manufactured in France,
- Power supply: 230 V,
- Emergency power supply: through 2 new built-in batteries: 12 V / 1.3 Ah,
- Frequency: 50/60 Hz,
- Power: 300 VA,
- Hill-Rom 900 is a multifunctional hospital bed designed to meet the needs of the entire medical team,
- Due to its functions, it is designed for use not only in general care, but also in post-operative care and the ICU,
- It facilitates the use of vital function monitoring devices and patient transport,
- Has buttons for selective function locking located on both sides of the bed,
- Infinite electric adjustment by wired remote control:
 - Bed height in the range: 40 - 78 mm (with mattress 57 - 95),
 - Back segment tilt angle in the range: 0 - 65º,
 - Trendelenburg and anti - Trendelenburg position: + 17 / - 17º,
 - Thigh segment tilt angle: 28º,
 - Foot segment tilt angle: - 3° to - 22°,
 - Chair position - EasyChair - comfortable and ergonomic sitting position,
 - Zero "0" position - the lowest and safest patient position,
 - Intelligent Auto Contour,
- Equipped with:
 - Safe, folding, long side rails with color-coded bumpers,
 - LINE-OF-SITE back segment angle indicator, and in Trendelenburg and anti-Trendelenburg positions,
 - Adjustable buckle covering the mattress (90 - 85 - 80 cm),
 - Accessory holders (drainage bags, urine bags, etc.), on both sides of the bed,
 - Levers for manual CPR on both sides,
 - Side control panels on both sides of the bed with locking/unlocking of individual functions,
 - Patient controls labeled with clear graphic symbols for intuitive operation of all its functions,
 - Removable head and leg tops equipped with blockades for stability during transport,
 - Integrated bed extension with a single movement allows to adjust the length within 20 cm (four positions: 5 cm, 10 cm, 15 cm, 20 cm) as required,
 - The bedding shelf allows you to conveniently place the linen when making the bed, preventing it from being deposited on the floor, improving work and hygiene (safe working load of 15 kg),
 - Intelligent night illumination (orange light) with lowest position indicator (green light),
 - 4 drip spaces in each corner of the bed,
 - 4 corner bumpers,
 - Brake lock alarm,
 - Central brake,
- X-RAY transparency,
- Mobile, on 4 wheels with Ø 15 cm, and additional self-steering 5th wheel,
- Includes:
 - Arjo Pentaflex anti-decubitus mattress, visible spots on the cover,
 - User manual in Polish (PDF),
- Dimensions:
 - Mattress: 87 x 197 x 17 cm,
 - Overall: 102 x 215 x 88.5/126 cm,
- Load capacity:
 - Safe: 185 kg,
 - Maximum: 220 kg,
- Weight: 153 kg,
- Cat 02,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.
WE SPEAK ENGLISH
HABLAMOS ESPANOL
WIR SPRECHEN DEUTSCH
ON PARLE FRANCAIS
WIJ SPREKEN DUTCH
Мы говорим по русски
Мы говорим по украински

